TECHNICAL SPECIFICATIONS FOR ANTISTATIC BUTADYL® GLOVE

| | ТҮРЕ | DESCRIPTION |
|-----------|-----------|-----------------|
| | MATERIAL | BUTADYL® |
| | COLOR | WHITE |
| | MODEL | EQ-AGB-Glove |
| | | HAND SPECIFIC |
| | HAND | OR |
| | STYLE | AMBIDEXTROUS |
| (FET LAND | THICKNESS | 30mm |
| | LENGTH | 32 INCHES |
| | BEAD | 0.25 INCH |
| | MARKING | MTI Corporation |

BUTADYL® CHARACTERISTICS

BUTADYL® IS A COST EFFECTIVE ALTERNATIVE TO HYPALON®, NEOPRENE OR BUTYL IN MANY APPLICATIONS. IT PERFORMS WELL IN SITUATIONS IN WHICH VAPORIZED HYDROGEN PEROXIDE (VHP) IS INVOLVED. IT ALSO OFFERS STRONG RESISTANCE TO PARTICULARLY HARMFUL AGENTS, SUCH AS: AROMATIC HYDROCARBONS, PETROLEUM SOLVENTS, AND LIPIDS. IT IS AN IDEAL MATERIAL THAT MAINTAINS CONSISTANCY WHEN FACED WITH PRESENCE OF OZONE OR SUN, AND IS NOT PRONE TO INDUCE ALLERGIC REACTIONS.

PERMEATION TEST RESULTS ON WHITE BUTADYL® GLOVES: ASTM F 739-99A

| ΤΕΣΤ CHEMICAL | BREAKTHROUGH PERMEATED MASS/UN | | STEADY STATE | |
|-------------------------------------|--|----------------------------|--------------------|--|
| | DETECTION TIME | AREA AT BREAKTHROUGH | PERMEATION RATE | |
| | MINUTES | (AVG.), UG/CM ² | (AVG.), UG/CM²/MIN | |
| ACETONE 99.5% | 58 MINUTES | 1097.9 | 814.0 | |
| ACETONITRILE 99.8% | 21 MINUTES | 514.3 | 2751.20 | |
| METHYL ALCOHOL (METHANOL) 99.84% | 90 MINUTES | 586.60 | 412.7 | |
| SULPHURIC ACID 37% | NO BREAKTHROUGH WAS DETECTED UP TO 240 MINUTES | -0- | -0- | |
| AMMONIUM HYDROXIDE 30% | NO BREAKTHROUGH WAS DETECTED UP TO 240 MINUTES | -0- | -0- | |
| HYDROFLUORIC ACID 51% | NO BREAKTHROUGH WAS DETECTED UP TO 240 MINUTES | -0- | -0- | |
| ANILINE 100% | 120 MINUTES | 123.8 | N/A | |

PHYSICAL PROPERTIES

| WHITE BUTADYL GLOVES | | | |
|------------------------|------|--|--|
| ULTIMATE ELONGATION, % | 417 | | |
| 100% MODULUS, PSI | 633 | | |
| 200% MODULUS, PSI | 1067 | | |
| 300% MODULUS, PSI | 1704 | | |
| TENSILE STRENGTH, PSI | 3848 | | |

GAS PERMIABILITY

| | GAS (REFERENCE) | AIR | HYDROGEN | NITROGEN | OXYGEN | HELIUM |
|----------|--------------------|-----|----------|----------|--------|--------|
| BUTADYL® | | 2.1 | 69.674 | 17.204 | 9.10 | 72.028 |

PERMEATION TEST RESULTS FOR CYTOTOXIC DRUGS, 4 HOUR DURATION

| | BREAKTHROUGH DETECTION TIME | STEADY STATE PERMEATION RATE | OTHER OBSERVATIONS |
|----------------------------|--------------------------------|---------------------------------|-----------------------|
| | Minutes | (Avg) Ug/cm²/min | |
| DACARBAZINE | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 10.0 MG/ML (10,000 PPM) | TO 240 MINUTES | | NO DEGRADATION |
| CARMUSTINE (BICNU) | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 3.3 MG/ML (3,300 PPM) | TO 240 MINUTES | | NO DEGRADATION |
| CYCLOPHOSPHAMIDE (CYTOXAN) | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 20.0 MG/ML (20,000 PPM) | TO 240 MINUTES | | NO DEGRADATION |
| DOXORUBICIN HYDROCHLORIDE | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 2.0 MG/ML (2,000 PPM) | TO 240 MINUTES | | NO DEGRADATION |
| 5-FLUOROURACIL | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 50.0 MG/ML (50,000 PPM) | TO 240 MINUTES | | NO DEGRADATION |
| CISPLATIN | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 1.0 MG/ML (1,000 PPM) | TO 240 MINUTES | | NO DEGRADATION |
| ETOPOSIDE (TOPOSAR) | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 20.0 MG/ML (20,000 PPM) | TO 240 MINUTES | | NO DEGRADATION |
| PACLITAXEL (TAXOL) | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 6.0 MG/ML (6,000 PPM) | TO 240 MINUTES | | NO DEGRADATION |
| ΤΗΙΟ-ΤΕΡΑ | NO BREATHROUGH UP | 0 | NO SWELLING, |
| 10.0 MG/ML (10,000 PPM) | TO 240 MINUTES | | NO DEGRADATION |

 ${}^{\rm Page}3$

PERMEATION TEST RESULTS FOR CYTOTOXIC DRUGS, 8 HOUR DURATION

| | BREAKTHROUGH DETECTION TIME | STEADY STATE PERMEATION RATE | OTHER OBSERVATIONS |
|----------------------------|--------------------------------|---------------------------------|-----------------------|
| | Minutes | (Avg) Ug/cm²/min | |
| CYCLOPHOSPHAMIDE | NO BREAKTHROUGH UP | 0 | NO SWELLING, |
| (CYTOXAN) 20.0 MG/ML | TO 480 MINUTES | | NO DEGRADATION |
| (20,000 PPM) | | | |
| DOXORUBICIN HYDROCHLORIDE | NO BREAKTHROUGH UP | 0 | NO SWELLING, |
| 2.0 MG/ML (2,000 PPM) | TO 480 MINUTES | | NO DEGRADATION |
| 5-FLUOROURACIL | NO BREAKTHROUGH UP | 0 | NO SWELLING, |
| 50.0 MG/ML (50,000 PPM) | TO 480 MINUTES | | NO DEGRADATION |
| METHOTREXATE | NO BREAKTHROUGH UP | 0 | NO SWELLING, |
| 25.0 MG/ML (25,000 PPM) | TO 480 MINUTES | | NO DEGRADATION |
| VINCRISTINE SULFATE | NO BREAKTHROUGH UP | 0 | NO SWELLING, |
| 1.0 MG/ML (2,000 PPM) | TO 480 MINUTES | | NO DEGRADATION |
| DAUNORUBICIN HYDROCHLORIDE | NO BREAKTHROUGH UP | 0 | NO SWELLING, |
| 5.0 MG/ML (5,000 PPM) | TO 480 MINUTES | | NO DEGRADATION |

PERMEATION RESULTS FOR VAPORIZED HYDROGEN PEROXIDE (VHP)

| CHEMICAL TESTED | TIME INTERVAL | AVERAGE BREAKTHROUGH DETECTION TIME | AVERAGE STEADY STATE PERM. RATE | COMMENTS |
|---------------------------|---------------|---|---|---|
| | | (SPECIMEN 1/2/3) Minutes | (SPECIMEN 1/2/3) (μg/cm²/minutes) | |
| | 1 HOUR | NO BREAKTHROUGH | N/A | NO SWELLING, DEGRADATION, OR SWELLING |
| HYDROGEN PEROXIDE, 35% | 5 HOURS | NO BREAKTHROUGH | N/A | NO SWELLING, DEGRADATION, OR SWELLING |
| | 10 HOURS | NO BREAKTHROUGH | N/A | NO SWELLING, DEGRADATION, OR SWELLING |
| | 20 HOURS | NO BREATHROUGH | N/A | NO SWELLING, DEGRADATION, OR SWELLING |

 $_{
m Page}4$

QUALITY CONTROL INSPECTION FLOW CHART



USAGE RECOMMENDATIONS

THIS GLOVE IS PARTICULARLY RECOMMENDED FOR:

- SODIUM HYDROXIDE 50%
- NITRIC ACID 10%
- POTASSIUM HYDROXIDE 50%
- PERCHLORIC AICD 60%

IT IS ALSO SUITABLE FOR HANDLING CHEMICAL PRODUCTS LIKE:

- ALCOHOLS
- MINERAL SPIRITS
- PENTANE, KEROSENE
- EDIBLE OILS AND ACIDS

THIS GLOVE IS NOT SUITABLE FOR HANDLING:

- KETONIC SOLVENTS (i.e. METHYL ETHYL KETONE)
- PLASTICISERS (i.e. PROPYLENE OXIDE)
- CHLORINATED SOLVENTS (i.e. METHYLENE CHLORIDE)
- AROMATIC SOLVENTS (i.e. XYLENE, STYRENE)
- PHENOL 90%

STORAGE AND MAINTENANCE

- BUTADYL GLOVES ARE PACKAGED FULLY EXTENDED, INDIVIDUALLY WRAPPED IN A POLYETHYLENE BAG.
- AFTER HANDLING, PLEASE REMOVE CHEMICAL RESIDUE FROM GLOVES WITH SOAPY WATER.
- BUTADYL GLOVES MUST BE RETURNED TO A POLYETHYLENE BAG AND STORED AT ROOM TEMPERATURE TO AVOID ALL RISKS OF DETERORATION AND DAMAGE.