





Solid-State Li-S R&D Equipment Line





Solid-State Li-S R&D Equipment Line

Dual Chamber Glovebox

Description

Dual Chamber Glove Box (94" x 31" x 36") with Gas Purification System (H₂O&O₂< 1ppm) and Openable Left Front Window –
[VGB-6-II-LD](#)



Dual Jars high Energy Ball Milling

Description

Dual Jars High Energy Ball Mill with Programmable Controller (Ar Gas Compatible) –
[MSK-SFM-3-II](#)



Muffle / Tube Furnace 2 in 1

Description

Compact Hybrid Furnace (4x4x4", 1000° C), Muffle/Tube 2-In-1 –
[KSL-1100X-S-H-LD](#)





Solid-State Li-S R&D Equipment Line

Compact Hydraulic Pellet Press

Description

15T Compact Hydraulic Pellet Press
for use in Glove-box - Made in
Europe –
[YLJ-15L](#)



Electric Roller Press

Description

4" Width Electric Roller Press with
Variable Speed - (Ar Glovebox
compatible) –
[MSK-HRP-MR100DC](#)




Mini Tap Casting / Doctor Blade Coating

Description

Mini Tape Casting Coater
(200Lx100W mm) with Heat-able
Vacuum Bed & Vacuum Pump at
100°C Max. –
[MSK-AFA-HC100](#)





Solid-State Li-S R&D Equipment Line

Hand-Held Disc Cutter

Description

Hand-Held Disc Cutter with One Set of Ring Dies with ID 8,15, 18, 20

mm –

[MSK-T-12](#)



Quick Assembly Split Coin Cell

Description

Quick Assembly Split Coin Cell for Battery Material testing, No

Crimper Needed –

[EQ-HSTC](#)

